

February 28th, 2019

To whom it may concern,

All channels are automatically calibrated to ensure linear and precise measurements as soon as the Limestone Technologies Inc. Paragon is plugged into the USB port of a computer. There are no calibration routines within the software except the Blood Pressure channel. A software calibration utility is incorporated into the software for this specific channel. Limestone Technologies recommends recalibrating the blood pressure channel if the instrument is being moved to different atmospheric pressures.



---

Tyler Buttle, Software Development and Support Manager, Limestone Technologies Inc.